

SJM Tailor® Flexible Ring + Band



[ease of implant for natural mitral and tricuspid repair]

SJM Tailor® Flexible Ring + Band

Flexible. Versatile. Designed to move with natural anatomy. Time-saving solutions for mitral and tricuspid valve repair.

The SJM Tailor® Flexible Ring + SJM Tailor® Flexible Band are fully flexible repair solutions designed to accommodate the natural movement of the annulus. This versatile system provides surgeons a wide range of options to meet any annuloplasty repair need – from a full ring that can be customized to any shape, to a standard pre-cut C-band.



SJM Tailor® Flexible Ring

- ✓ Provides anterior and posterior support while allowing dynamic annular motion.
- ✓ Full, flexible ring may be customized to any shape addressing specific patient needs.


SJM Tailor® Flexible Band

- ✓ Provides posterior support while allowing dynamic annular motion.
- ✓ Convenient pre-cut C-band for time-saving simplicity.

Flexibility in just one step.

Delivering ease of implant.

In addition to the customization and convenience of the SJM Tailor® Flexible Ring + Band, the St. Jude Medical delivery system employs a one-step process designed to save time and reduce complexity during surgery. Key implantability benefits are achieved through product features such as:




- Large target area eases suturing.
- Open ring holder facilitates visualization and maintains band shape while suturing.
- One-step push button release simplifies the implantation process.
- Secure ergonomic handle quickly attaches and detaches to save time.
- Malleable sizers and extension handle bend easily to aid maneuverability.



Extension handle can be quickly and securely attached for added length during minimally invasive procedures.


Inspired by the flexible motion of the annulus.

The SJM Tailor® Flexible Ring + Band provide the options needed for both mitral and tricuspid valve repair. With a fully flexible design that allows dynamic annular motion, the customizable ring and convenient C-band offer simple, easy-to-use alternatives for annuloplasty repair. As part of the broad St. Jude Medical Repair Solutions portfolio of products, the SJM Tailor® Flexible Ring + Band deliver the quality design surgeons expect from the worldwide leader in cardiac surgery solutions.



flexible.versatile.

SJM Tailor® Flexible Ring + Band



SJM Tailor® Flexible Ring + Band

Reference Dimensions

Ring Model Number	Band Model Number	Size	Inside Dimension	Outside Dimension	Inter-trigonal Dimension
TARP-25	TAB-25	25 mm	29 mm	34 mm	25 mm
TARP-27	TAB-27	27 mm	31 mm	37 mm	27 mm
TARP-29	TAB-29	29 mm	34 mm	40 mm	29 mm
TARP-31	TAB-31	31 mm	36 mm	42 mm	31 mm
TARP-33	TAB-33	33 mm	39 mm	45 mm	33 mm
TARP-35	TAB-35	35 mm	41 mm	46 mm	35 mm

Accessories

Model Number	Contents
TAR-505	SJM Tailor Annuloplasty Ring Sizer Set
	One (1) malleable holder handle
	One (1) extension handle
	Six (6) sizers (25, 27, 29, 31, 33, 35)
	One (1) autoclavable tray for storage of components
	One (1) autoclavable tray cover
HH-05	Replacement Holder Handle
EX-05	Replacement Extension Handle

VISIT OUR WEB SITE AT www.sjm.com

CAUTION: FEDERAL LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.
Please see the physician's manual for a full description of indications, contraindications, side effects, precautions, warnings and instructions for use.

Corporate Headquarters St. Jude Medical, Inc., One Lillehei Plaza, St. Paul, Minnesota 55117 USA **24-Hour Technical/Professional Consultation** (800) 328-9634 (USA) (651) 483-2000 Fax: (651) 482-8318 **Customer Service** (800) 544-1664 (USA) (651) 490-4410 Fax: (651) 481-7702 **European Headquarters** St. Jude Medical Europe, Inc., The Corporate Village, Avenue Da Vinci laan, 11-Box F1, B-1935-Zaventem, Belgium **Customer Service** Tel: 32-2-774-68-11 Fax: 32-2-772-83-84 **Asian Headquarters** St. Jude Medical Hong Kong Limited, Room 2705-2708, China Merchants Tower, Shun Tak Centre, 168-200 Connaught Road, Central, Hong Kong Tel: (852) 2996-7688 Fax: (852) 2956-0622.

SJM Tailor is a trademark of St. Jude Medical, Inc. US Pat # 6,174,332; patent and foreign patents pending. St. Jude Medical Cardiac Surgery Division ©2006 St. Jude Medical, Inc. All rights reserved.
Printed in the USA. ITEM1896/0206/5M/EN/YC

QUALITY. RELIABILITY. INTEGRITY.

 **ST. JUDE MEDICAL**