

Prodigy MRI™ Patient Programmer with BurstDR™ Stimulation*

Model 3856

Compatibility

- Prodigy MRI™ IPG Model 3772

Hardware Specifications

Specifications	Metric	Imperial
Height	6.80 cm	2.7 in
Length	10.77 cm	4.2 in
Thickness	2.60 cm	1.0 in
Weight	128.0 g	4.6 oz
Batteries	3 AAA	

Accessories

Description	Model #
Patient magnet	1210
Communication wand	1232
Battery pack	1254
Programmer case	1272



*BurstDR™ stimulation, patented technology exclusively from St. Jude Medical, is also referred to as Burst stimulation in clinical literature.

CE marking does not necessarily indicate regulatory approval status for all markets. Please refer to the instructions for use for a full listing of indications, contraindications, warnings and precautions.

Rx Only

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Indications for Use: Spinal cord stimulation as an aid in the management of chronic, intractable pain of the trunk and limbs. **Contraindications:** Patients who are unable to operate the system or who have failed to receive effective pain relief during trial stimulation. **Warnings/Precautions:** Diathermy therapy, implanted cardiac systems, magnetic resonance imaging (MRI), explosive or flammable gases, theft detectors and metal screening devices, lead movement, operation of machinery and equipment, postural changes, pediatric use, pregnancy, and case damage. Patients who are poor surgical risks, with multiple illnesses, or with active general infections should not be implanted. **Adverse Effects:** Painful stimulation, loss of pain relief, surgical risks (e.g., paralysis). Clinicians manual must be reviewed for detailed disclosure.

Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL and the nine-squares symbol are trademarks and service marks of St. Jude Medical, Inc. and its related companies. © 2016 St. Jude Medical, Inc. All Rights Reserved.

St. Jude Medical Inc.

Global Headquarters
One St. Jude Medical Drive
St. Paul, MN 55117
USA
T +1 651 756 2000 | F +1 651 756 3301

St. Jude Medical S.C., Inc.

Americas Division
6300 Bee Cave Road
Bldg. Two, Suite 100
Austin, TX 78746
USA
T +1 512 286 4000 | F +1 512 732 2418

SJM Coordination Center BVBA

The Corporate Village
Da Vincilaan 11-Box F1
B-1935 Zaventem, Belgium
T +32 2 774 68 11 | F +32 2 772 83 84

St. Jude Medical Brasil Ltda.

Rua Itapeva, 538
5º ao 8º andares
01332-000 – São Paulo – SP
Brazil
T +55 11 5080 5400 | F +55 11 5080 5423

St. Jude Medical (Hong Kong) Ltd.

Suite 1608, 16/F Exchange Tower
33 Wang Chiu Road
Kowloon Bay, Kowloon
Hong Kong SAR
T +852 2996 7688 | F +852 2956 0622

St. Jude Medical Japan Co., Ltd.

Shiodome City Center 15F
1-5-2 Higashi Shinbashi, Minato-ku
Tokyo 105-7115
Japan
T +81 3 6255 6370 | F +81 3 6255 6371

St. Jude Medical Australia Pty, Ltd.

17 Orion Road
Lane Cove, NSW 2066
Australia
T +61 2 9936 1200 | F +61 2 9936 1222

SJMprofessional.com